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The Maltese cross tube is used to demonstrate the linear motion of electrons. When the tube is set up in the magnetic field of a pair of Helmholtz coils (555 06), it becomes possible to conduct experiments on electron optics.

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1 Safety notes

- Heating voltages greater than 6.3 V AC reduce the service life of the tube. Voltages greater than 7 V AC can destroy the tungsten filament!
- Do not jar or jerk the tube while the filament is glowing or directly after switching off the heating (danger of breaking the hot filament!)
- Never connect more than one connector each to the 4-mm plug pins. Two plugs connected on top of each other may cause the glue joint of the plastic cap to break.

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2 Description (see Fig. 1), technical data

- ① Electron gun, consisting of directly heated cathode (tungsten filament) and cylindrical anode
- ② 4-mm socket pair for connecting the hot-cathode heating voltage
- ③ 4-mm plug pin connected with anode
- ④ Maltese cross, aluminum, as shadow body
- ⑤ Luminescent screen
- ⑥ 4-mm plug pin connected with Maltese cross

The electrode system ① of the tube consists of a tungsten wire "hairpin" cathode and a cylindrical anode. The two terminals of the tungsten cathode are connected to 4-mm sockets ② in a plastic cap at the tube neck. The anode terminal is led out via a 4-mm plug pin ③ on the side of the tube neck.

The cathode generates a divergent electron beam which is incident on the luminescent screen ⑤ attached on the opposite wall of the glass tube. An aluminum Maltese cross ④ is located

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
between the electrode system and the luminescent screen; the cross is connected to plug pin ⑥ on the top plastic cap.

Technical data

Anode voltage:	2 kV DC to max. 5 kV DC
Heating voltage:	approx. 6 V AC ... 6.4 V AC
Recommended:	6.3 V
Diameter of glass tube:	approx. 13 cm
Tube length overall:	approx. 25 cm

3 Operating equipment

- Stand for electron tubes (555 05)
 Voltage supply for the tube
- Anode and heating voltage from high-voltage power supply (521 70) or
 - Anode voltage from high-voltage power supply 10 kV (522 37)*
 Heating voltage (6 V AC) e.g. from transformer 6 V; 12 V; 30 W (522 73)

4 Experiment examples

Fig. 2

4.1 Linear propagation of cathode rays

Set up the circuit shown in Fig. 2. The anode of the beam generation system and the metal cross must have the same potential. First, turn on only the heating voltage: we observe that the visible light emitted by the hot cathode projects a sharp shadow of the cross on the luminescent screen.

Now switch on the anode voltage to generate a divergent electron beam. You can now observe that a sharp shadow caused by the charged particles appears on the screen, fully congruent with the first shadow. Adjust the voltage in the beam generation system to obtain a suitable intensity of the luminescence.

This experiment demonstrates that the charges (cathode rays) produce a shadow in exactly the same way as visible light, and must thus propagate linearly.

4.2 Deflection in a magnetic field

While the tube is in operation, place one or more magnets (e.g. 510 48) in the vicinity of the tube and observe the following phenomena:

1. The shadow is shifted.
2. The extent of the shift depends on both the magnetic field strength and the acceleration voltage which is applied to the beam generation system.
3. The shadow is not deflected in the direction of the magnetic field.

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
If we relate the direction of deflection, the direction of the field and the direction of motion of the charges using the right-hand rule, we find that:

Cathode rays seem to behave in a magnetic field analogously to currents in conductors.

4.3 Effects of electrostatic charges

Disconnect the metal cross from the anode potential and observe how the shadow is distorted. A negative charge accumulates on the cross. When the equilibrium state is reached, this inhibits the accumulation of additional charges. Cathode rays which pass close to this opposing field are deflected. This distorts the shadow as shown in Fig. 3. When the cross is connected to the cathode potential, the image is distorted so much that it extends beyond the edges of the luminescent screen.

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4.4 Introduction to electron optics

Mount the tube in the retaining straps of the stand so that the screen is surrounded by one of the two Helmholtz coils (555 05). In contrast to its normal position, the Maltese cross tube is mounted in the clamp jaws so that in this experiment the base of the tube extends through the other Helmholtz coil.

Increasing the magnetic field, i.e. the coil current, causes the shadow to rotate, contract to a small spot, then turn further and grow in the opposite direction. Changing the anode voltage makes possible further changes in the shadow image.

Thus, cathode rays and deflection fields can be used to enlarge electronic screen images in a manner analogous to an optical lens system.